## English Original

### **Device Description:**

The Parodi Anti-Emboli System (PAES) protects the brain from emboli released during interventional procedures by reversing the direction of blood flow in the internal carotid artery. It is deployed proximal to the lesion being treated and either passive or active flow reversal maneuvers may be used. The PAES consists of the 10 F OD Parodi Anti-Emboli Catheter™ (PAEC™), the Parodi External Balloon™ (PEB™) and the Parodi Blood Return Set™ (PBRS™). The devices are supplied with two 3 ml syringes and two 4-way stopcocks for balloon inflation, and a peel-away introducer sheath for PEB insertion.

Blood flow is reversed when the PAEC and PEB occlusion balloons are inflated, isolating the internal carotid artery and creating a pressure gradient through the catheter. The proximal exit port of the PAEC is connected to the PBRS, which returns filtered arterial blood to the femoral vein. When needed, flow reversal can be actively supplemented by aspiration through the guiding catheter using the 4-way stopcock on the PBRS.

At the completion of the procedure the balloons are deflated, restoring antegrade blood flow, and the PAEC is removed.

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